



Carum

True Hospital Pump

The Carum hospital-grade breast pump has many benefits: it is powerful while also being ideally suited for use in sensitive situations. It is also unbeatable in the double pump sector thanks to its dual-piston technology.



INDIVIDUAL PHASES

Carum is the first hospital pump to allow the cycle and vacuum intensity settings to be controlled independently from one another in both stimulation and expression mode.



SENSITIVE PROGRAMME

In the "Sensitive Programme" the Carum operates with a barely detectable vacuum and cycle, thus enabling gentle pumping if nipples are painful and sore. The settings can also be adjusted to the sensitivity of the mother in small steps.



CLOSED SYSTEM

The proven "Vacuum Seal" technology serves as a 100 % barrier against bacteria and ensures that no breast milk or impurities can enter the connecting tube or pump. This provides reliable protection against contamination of the breast milk and thus prevents infection.

integrated
double bottle
holder



closed
system



strong suction,
250 mmHg single or
double pumping



integrated battery -
4 hour battery life



timer to keep
track of your
pump session



dimming
function



CARUM

WELCOME

ARDO

smooth surface
is easy to clean



suction & cycle
independently
adjustable,
64 settings



Sensitive Programme -
low level vacuum for
sore nipples



let-down button for
2-phase function,
similar to the sucking
behavior of baby



(844) 411-ARDO



sales@ardo-usa.com



www.ardo-usa.com